


For further information:

P: 06 878 8109 ext. 5761

W: www.tewhatauora.govt.nz

Health New Zealand
Te Whatu Ora

January 2025




Shower Stool

Instructions for setting up
a shower stool at home

Health New Zealand
Te Whatu Ora

Instructions for setting up a shower stool at home

To adjust height:



- Push in the metal buttons on the leg.
- Set legs at the height to suit the person using the shower stool.
- Click the buttons into the correct holes. Make sure buttons are completely popped out.
- Make sure all four legs are the same height.
- One leg has extra holes - this is used to adjust the stool for uneven surfaces.

To use stool:



- Make sure feet are flat on floor and hips are above the knees.
- Always use two hands to push off the frames arms as it can tip.
- Stool can be used in a cubicle or shower, or for sponge washing at the basin.
- Stool can be used with non-slip bath mat if surfaces are slippery.
- **MAXIMUM WEIGHT CAPACITY IS 120KG.**
DO NOT USE STOOL IN THE BATH OR SHUB

Additional information from your Occupational Therapist:

If you have any questions or require further advice please ring the Hawke's Bay Hospital Equipment Loan Store
P: 06 878 8109 ext. 5761